

## **Solution Phase Combinatorial Chemistry in Drug Discovery**

**Ischia Advanced School of Organic Chemistry  
September 29-October 4 1996**

**Richard Storer, GlaxoWellcome Research and Development**

GlaxoWellcome

## **Industry Requirements**

- **First and best**
- **Reduced timescales**
- **More targets to evaluate**


Development

GlaxoWellcome

## Compound Libraries in Drug Discovery

- Primary libraries for hit generation
  - pools or discretés
- Hits to Leads
- Lead optimisation
  - discretés

GlaxoWellcome



## Solution Phase Chemistry

- Wide range of established chemistry
- Short lead-in time
- Ideal for short sequences
- No cleavage needed
- Direct link with screens
- High throughput
- Quantity

GlaxoWellcome

## Library Size

Dimer ----- Trimer ----- Tetramer

increasing numbers  
increasing diversity

but: increasing steps  
increasing M.Wt.

Automated

S<sub>x</sub>  
C<sub>y</sub>-B<sub>z</sub>

GlaxoWellcome

## Amide Trimer Synthesis

Carboxylic acid-----Aminoacid-----Amine

80

50


40

Amide bond chemistry  
 160,000 compounds  
 Pool size of 40  
 Quantity: 25 micromoles  
 Semi-automated synthesis


GlaxoWellcome

## Manual/Automated Approach

- Manual preparation of feedstock dimers
- Automated synthesis of final trimers



- Decoding - reverse direction synthesis
- Maintain automated synthesis of final trimers




## Chemistry Development

- Discrete synthesis
  - all steps, full analysis
  - monomer selection
- Pooled synthesis
  - nmr, hplc, ms, lc-ms

GlaxoWellcome

## Screening Controls

- Standards using automated procedure
- Inclusion of known lead / active compounds



- By-product screening

Robotic Chemistry


GlaxoWellcome

## First Generation Robotics

- Solution dispensing
- Multiple reaction vessels
- Variable reagent addition
  - order
  - volume
- Heating capacity
- Solvent removal

GlaxoWellcome


## Library Synthesis



GlaxoWellcome


## Synthetic plate design

- Direct parallel to screening mtp - 80 well plate
- Single aminoacid component on plate
- 80 acid array on each plate
- 40 amine components in each well
- Compound identity by plate/well address




GlaxoWellcome

## Decode Synthesis




GlaxoWellcome

## HIV Proteinase Inhibitor



IC<sub>50</sub> = 130 µg/ml


Sample decoded to individual compounds



GR204887X IC<sub>50</sub> = 10 µg/ml


GlaxoWellcome

## NK1 Antagonist



IC<sub>50</sub> = 8 µg/ml

Sample decoded to individual compounds



GR 217869X pKi 7.3

AX file

GlaxoWellcome




## Profiling of libraries

- Descriptor based
  - Charges, H-bonds, size, rings, flexibility
- Applications
  - Diversity within libraries
  - Comparison between libraries
  - "Hole filling" in collections
  - Target profiling

GlaxoWellcome

## Library comparison with drug files




GlaxoWellcome

## Conclusions


- Highly successful exercise
- Leads generated in many areas
- Projects established
- Long-term screening resource generated
- General trend information without decoding
- Decoding causes problems
  - success rate
  - timescales

GlaxoWellcome


## Future Directions




Total compounds: 400K




58K



182K



96K



38K

GlaxoWellcome

## Future Directions

- **Extend range of chemistries**
  - dimer formation / trimer formation
- **Discrete synthesis / pooling option**
- **Extend robotic capability**

GlaxoWellcome

## Chemistry Requirements

- **Clean and high yielding**
  - one major pathway
  - few (no) by-products
- **Versatile and robust**
- **Combinatorial ( A + B  $\longrightarrow$  C )**
- **Mild conditions**
- **Availability of building blocks**
- **Suitable products**

GlaxoWellcome

## Other Developed Solution Chemistry

- 1,4 Additions
- Aldols
- Cycloadditions
- Heteroaryl chloride displacement
- 5,6,7 ring heterocycles
- Hydrazones
- N-alkylation
- Oximes
- Reductive amination
- S-alkylation

GlaxoWellcome



## Heterocycle Libraries

- Chemically efficient
- Provide good templates
- Ubiquitous in drug molecules

Benzimidazoles	Pyrazolones
Benzodiazepines	Pyrimidines
Imidazopyridines	Thiazoles
Oxadiazoles	Triazoles


GlaxoWellcome

## Oxadiazole Libraries



GlaxoWellcome


## Thiazole Dimer Libraries



1440 discrete compounds in each series

S.P. Watson et al., *Bioorg Med Chem Lett*, 1996, 6, 1406


## Three Component Libraries



- Large libraries from modest monomer numbers
- 11,520 products from 104 monomers
- Lipophilic, H-bonding, acidic and basic substituents


GlaxoWellcome

## Pyrazolone Condensation Trimers



GlaxoWellcome


## Anhydride Trimers



GlaxoWellcome


## Bradykinin Diamide Hit

- One hit from 4000 mixture library
- Decoded to single active component



GR213548X pKI = 5.8


## Hit to Lead Focused Library



GlaxoWellcome

## Histamine H3

- Project Lead
- Potential metabolic problem - imidazole




GR175737

J.W.Clitherow et al., *Biorg Med Chem Lett*, 1996, 6, 833

GlaxoWellcome




## Project Route



GlaxoWellcome


## Initial Robotic Route



- Reaction failed when tertiary amines in R group

GlaxoWellcome

## Modified Robotic Route




- Extended to discrete library of 1600 samples

GlaxoWellcome

## Influenza Sialidase

- GG167 - Treatment for influenza A and B
- Follow up - simplified structures



PW Smith et al., *Biorg Med Chem Lett*, 1996, 6, 1805

GlaxoWellcome

## **New Generation Robotics**

- **Microtitre blocks**
- **Agitation**
- **Heating**
  - conventional oven
  - microwave oven
- **Vacuum centrifugation**
- **Solution phase purification**
  - aqueous extraction
  - filtration

GlaxoWellcome

## **Future Directions**

- **Low-tech equipment for modest numbers**
- **Equipment for each laboratory**
- **Array synthesis of larger (gm) quantities**
- **Parallel purification**
- **Solid phase arrays in lead optimisation**
- **Supported reagents**
- **Effective use of the technologies**

GlaxoWellcome

## Project Considerations

- Wider choice of leads available
- Generate leads at the appropriate time
- Define preferred series
- Screen capacity must parallel chemical capacity
- Arrays for optimisation and evaluation
- Capacity to make simultaneous changes
- Results on x80 give better strategic view
- Less justification for individual compounds

GlaxoWellcome

## Future Issues

- Solution or solid phase
  - as appropriate in a particular case
- Libraries
  - smarter not bigger
  - expandable chemistry
- Monomer availability
- Process chemistry

GlaxoWellcome

## Acknowledgements

Nick Bailey  
Mark Bamford  
John Bradshaw  
Paul Bronnert  
Darko Butina  
Stel Christodoulou  
Tony Cooper  
Martyn Deal  
Tony Dean  
Andrea Gore  
Duncan Judd  
Chris Marshall  
Judith Menzies  
Andy Merritt  
David Middlemiss  
Mark Ness

Chris Newall  
Tony Parkhouse  
Nigel Parr  
Neil Pegg  
Adrian Pipe  
Jan Scicinski  
Peter Seale  
James Smart  
Steve Spooner  
Jun Uchida  
Peter Ward  
Steve Watson  
Hazel Weldon  
Julia Widdowson  
Richard Wilson  
Chris Wooley

GlaxoWellcome