

Designing Degradable Implant Materials

Peter Neuenschwander and Ulrich W. Suter
Department of Materials, Institute of Polymers, ETH Zürich

with the invaluable assistance and support of

Gianluca Ciardelli
Martin Colussi
Marco Camus
Carmelina Grob-Pisano
Thomas D. Hirt

Olivier M. Keiser
Phillippe Kern
Kazushige Kojima
Andreas Lendlein

Sandro Matter
Bashar Saad
Paul Schmutz
Rémy C. Stoll
Kimberly Willis



Collaborating Medical and Veterinary Units

Patrick Aebischer

Centre Hospitalier Universitaire Vaudois &
Ecole polytechnique fédérale de Lausanne

Günter Burg

Dermatology, University of Zürich

Hans Gerber

AO Center, Association for the Study of
Internal Fixation, Davos

Brigitte von Rechenberg

Veterinary Surgery, University of Zürich

Georg K. Uhlschmid

Surgery, University of Zürich

Gregor Zünd

Cardiovascular Surgery, University of Zürich



Why New Medical Implant Materials? What Should Be Their Characteristics?

- most of the polymers currently used in implants were not initially designed for medical use —
rapidly increasing awareness of potential risks
- each implant application calls for a specific set of properties; hence, **materials are needed in which the most important characteristics can be modified independently from the others (materials with orthogonal properties)**



Blockcopolymers: Materials with Orthogonal Properties


Natta, Mazzanti, Crespi & Moraglio *La Chimica e l'Industria* **1957**, 39, 275
"Stereoblock copolymers of polypropylene (PP)"




Materials Properties — Findings:

- The isotactic blocks crystallize, the atactic blocks form a continuous amorphous phase (T_g as in atactic PP)
- The crystal T_m depends on the length of the isotactic blocks
- Between T_g and T_m , the material is a (thermoplastic) **ELASTOMER**
- The material's modulus is determined by the number of crystals (often by the volume fraction of crystalline blocks)

**crystalline α -form
of isotactic PP**



Approach to Degradable Block-Copolymers



hard segments
form domains
that are crystalline
or glassy

(→ THERMO-MECHANICS)

soft segments
form domains
that are amorphous
and "rubbery"

(→ DEGRADE FIRST)

● junction units

macrodiol 1 + macrodiol 2 + junction unit

PHB-diol

diisocyanate or
diaciddichloride

block-copolymer


copolyesterurethane or
copolyester

Polymer

ETH

Materials

Hard Segments: Telechelic PHB-diols by Transesterification of Bacterial PHB with a Diol



catalyst
dibutyl-Sn(II)-dilaurate

Polymer

ETH


Materials

T. Hirt


PHB-diol (M_n 2'500)
+ Poly(ϵ -caprolactone) (M_n 2'500)
+ (S)-2,6-Diisocyanato methylcaproate

PHB-diol (hard segment) [wt-%]	PCL-diol (soft segment) [wt-%]	tensile modulus [MPa]	tensile strength [MPa]	elongation at break [%]	T_m [°C]
29	58	60	6	830	116
44	44	200	7	610	122
69	21	500	11	21	134
92	0	1260	14	2	131

**The Modulus is Mainly Dependent on
 the Fraction of Hard Segment**
 hard segment: PHB (M_n 2'500)



Control of Degradation Rates: Introduction of "Weak Links"



hard segments
(→ THERMO-MECHANICS)

soft segments
(→ ELASTOMERIC
PROPERTIES)

"weak links" (readily
hydrolyzable groups)
(→ DEGRADATION RATE)

"weakest link":
glycolyl glycolate




Polymer




A. Lendlein

Soft Segments by Ring-Opening Polymerization from Diglycolide and ϵ -Caprolactone



A. Lendlein


Hydrolysis of Block-Copolyesterurethanes with 44 wt-% PHB and Different Soft Segments



pH 7 — 37 °C

Fluorescent Labelling of Telechelic PHB-Diol with a Dye

The dansyl-labelled blocks form hard-segment crystals that can be traced by their fluorescence




Phagocytosis of Fluorescent-labelled PHB Particles in Rat Macrophages

macrophage cell line J774 — 2 $\mu\text{g/ml}$ DPHB particles
optical micrographs after 4 hours

phase-contrast image

fluorescence image



polymer

ETH


Materials

B. Saad

Macrophage Cell Cultures that Were Exposed to PHB Particles as They Result from Significantly Degraded DegraPol[®]



PHB (2 μg PHB/ml)



PHB

primary labelled
peritoneal
macrophages

polymer
M. Müller

B. Saad, G.K. Uhlschmid

polymer

ETH

Materials

Phagocytosis of Fluorescent-labelled PHB Particles in Co-cultures of Hepatocytes and Kupffer Cells

1 $\mu\text{g}/\text{ml}$ of DPHB particles, after 4 hours

phase-contrast image

fluorescence image



hepatocytes


— 10 μm

Kupffer cells phagocytose PHB particles, hepatocytes are unaffected

Journal of
MEM
Materials

P. Schmutz



Functionalization of Copolyesterurethanes: Adding Carboxylic Side Groups



COOH



anchor unit




10,11-dihydroxy-undecanoic acid

Journal of
MEM
Materials

G. Ciardelli

Reaction of Polymer-bound COOH with "Any" Amine Through Activation with *p*-Nitrophenol

1. Formation and isolation of the
2. Reaction with amines — no further activation necessary




© 1999 EMM Materials, Inc.

EMM
Materials

G. Ciardelli


Nerve Conduits — Operating Principle



© 1999 EMM
Materials

P. Aebischer

Nerve Conduits Made from DegraPol® with Varying PHB Content



ETH
Materials

R. Stoll

Nerve Conduits Made from DegraPol



Whole sciatic nerve fully regenerated in a DegraPol nerve conduit after 12 weeks.



Histological section of a nerve regenerated neural tissue in a DegraPol nerve conduit with myelinated axons grouped into nerve fascicles.

Number of axons is not significantly different than in a normal sciatic nerve

ETH
Materials

P. Aebischer

Production of Open-pored Scaffolds

DegraPol
solution in
1,4-Dioxane



phase segregation into
polymer-rich and
polymer-poor phase

3D structure of
polymer-rich
phase is fixed




open-pored
DegraPol scaffold




DEGRA POL
Materials

S. Matter

Porous Foams with Different Pore Size



"normal"




the smallest pores
practical

DEGRA POL
Materials

S. Matter, M. Camus

Porous Foams with Different Pore Sizes in One Object



structures with small and large pores can be joined through an abrupt transition region, without "skin" between the regions


ETH
materials

M. Camus

Porous Tubes from DegraPol[®] (Obtained by Extrusion of a Gel in 1,4-Dioxane into Liquid N₂ and Freeze Drying)



optical micrograph
(tube appears white because of light scattering)




SEM

Dr. M. Müller
ETH ZH 185516/87/ETH

ETH
materials

R. Stoll, Ph. Kern

Application: Foams as Tissue-culture Substrates



DegraPol® Scaffold disks, diameter 10 mm, thickness 0.5 or 1 mm
24 pieces packed in a 24-well tissue-culture plate

ETH
Materials

M. Camus

An Application: Cartilage Repair




adult pig chondrocytes, cultured on DegraPol foam for 2 weeks, then
stained with Alcian blue for Collagen II

→ the correct intercellular matrix is produced

ETH
Materials

B. von Rechenberg, B. Saad

An Application: Cartilage Repair



adult pig chondrocytes, cultured on Degra-Pol foam for 4 weeks, then stained with toluidine blue to visualize cells
 → **chondrocytes proliferate in the interstitial spaces of the foam**

DEGRA-POL
 Materials

B. von Rechenberg, B. Saad

An Application: Cartilage Repair



plug of Degra-Pol sponge, seeded with chondrocytes and cultured for 2 weeks, then glued to a cortical plug of hard bone with liquid bone cement.

The plug will be inserted in a cortical hole drilled into the bone and be fixed with liquid bone cement.

DEGRA-POL
 Materials

B. von Rechenberg

An Application: Cartilage Repair



our e-actant plug with a 100-
gig chondrocyte-culture.
DegraPol sponge (poly-
degradable variety),
implanted in the transected
bovine femoral (line) of a pig
for 3 weeks

The implants form a smooth surface with the joint, no abrasion,
also at the counterpart (femour condyles). DegraPol foam is yellow
because of small contamination with blood

DEGRA
Materials

B. von Rechenberg

Conclusions

- Multi-Blockcopolymers allow the design of a large variety of polymeric, tissue-compatible materials with a number of continuously tunable, orthogonal properties
- Independently variable are: softening point, mechanical moduli, degradation rates, functionalizability
- Functionalization can be achieved through attachment of any amine-terminated compound
- The polymeric materials described here (DegraPol[®]) have been extensively tested in vitro and in vivo and are highly tissue compatible
- The DegraPol[®]-foam structures are the first available elastomeric tissue-culture scaffolds

DEGRA
Materials